

Product datasheet for **GA206270**

Mouse Trem3 activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	Trem3 CRISPRa kit - CRISPR gene activation of mouse triggering receptor expressed on myeloid cells 3
Symbol:	Trem3
Components:	GA206270G1, Trem3 gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074) GA206270G2, Trem3 gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074) GA206270G3, Trem3 gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074) 1 CRISPRa-Enhancer vector, SKU GE100056 1 CRISPRa scramble vector, SKU GE100077
UniProt ID:	Q9JKE1
Synonyms:	BB134760
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online »](#)