

## Product datasheet for **GA206279**

### Mouse Shank3 activation kit by CRISPRa

#### Product data:

Product Type:	CRISPRa Kits
Description:	Shank3 CRISPRa kit - CRISPR gene activation of mouse SH3 and multiple ankyrin repeat domains 3
Symbol:	Shank3
UniProt ID:	<u><a href="#">Q4ACU6</a></u>
Synonyms:	AI841104; Shank3b; SHANK3c-3; SHANK3c-4
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

#### Product images:



[View online >](#)