


## Product datasheet for **GA206288**

### Mouse Gpr180 activation kit by CRISPRa

#### Product data:

|               |                                                                                                                                                                                                                                                                                                               |
|---------------|---------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| Product Type: | CRISPRa Kits                                                                                                                                                                                                                                                                                                  |
| Description:  | Gpr180 CRISPRa kit - CRISPR gene activation of mouse G protein-coupled receptor 180                                                                                                                                                                                                                           |
| Symbol:       | Gpr180                                                                                                                                                                                                                                                                                                        |
| Components:   | GA206288G1, Gpr180 gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074)<br>GA206288G2, Gpr180 gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074)<br>GA206288G3, Gpr180 gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074)<br>1 CRISPRa-Enhancer vector, SKU GE100056<br>1 CRISPRa scramble vector, SKU GE100077 |
| UniProt ID:   | <u><a href="#">Q8BPS4</a></u>                                                                                                                                                                                                                                                                                 |
| Synonyms:     | AB041545; E130016I23Rik; Itr                                                                                                                                                                                                                                                                                  |
| Format:       | 3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)                                                                                                                                                                                                                                       |

#### Product images:



[View online »](#)