

Product datasheet for **GA206294**

Mouse Eid1 activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	Eid1 CRISPRa kit - CRISPR gene activation of mouse EP300 interacting inhibitor of differentiation 1
Symbol:	Eid1
Components:	GA206294G1, Eid1 gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074) GA206294G2, Eid1 gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074) GA206294G3, Eid1 gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074) 1 CRISPRa-Enhancer vector, SKU GE100056 1 CRISPRa scramble vector, SKU GE100077
UniProt ID:	<u>Q9DCR4</u>
Synonyms:	2610002K20Rik; Cri1; EID-1; ORF12; PNAS-22; PTD014
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online »](#)