

Product datasheet for **GA206298**

Mouse Trpm7 activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	Trpm7 CRISPRa kit - CRISPR gene activation of mouse transient receptor potential cation channel, subfamily M, member 7
Symbol:	Trpm7
UniProt ID:	Q923J1
Synonyms:	2310022G15Rik; 4833414K03Rik; 5033407O22Rik; CHAK; CHAK1; Ltrp7; LTrpC-7; Ltrpc7
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online >](#)