

## Product datasheet for **GA206307**

### Mouse *Tssk3* activation kit by CRISPRa

#### Product data:

Product Type:	CRISPRa Kits
Description:	Tssk3 CRISPRa kit - CRISPR gene activation of mouse testis-specific serine kinase 3
Symbol:	Tssk3
UniProt ID:	<a href="#">Q9D2E1</a>
Synonyms:	1700014N07Rik; 4930594I21Rik; Stk22c; Stk22d; Tsk-3; Tssk-3; Tssk3b
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

#### Product images:



[View online >](#)