


Product datasheet for **GA206308**

Mouse Tdh activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	Tdh CRISPRa kit - CRISPR gene activation of mouse L-threonine dehydrogenase
Symbol:	Tdh
Components:	GA206308G1, Tdh gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074) GA206308G2, Tdh gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074) GA206308G3, Tdh gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074) 1 CRISPRa-Enhancer vector, SKU GE100056 1 CRISPRa scramble vector, SKU GE100077
UniProt ID:	<u>Q8K3F7</u>
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online »](#)