

Product datasheet for GA206320

Mouse Smpd3 activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	Smpd3 CRISPRa kit - CRISPR gene activation of mouse sphingomyelin phosphodiesterase 3, neutral
Symbol:	Smpd3
UniProt ID:	Q9JLY3
Synonyms:	4631433G07Rik; AI427456; AW537966; fro; nSMase2
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online >](#)