


## Product datasheet for **GA206330**

### Mouse Sqor activation kit by CRISPRa

#### Product data:

|               |                                                                                                                                                                                                                                                                                                         |
|---------------|---------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| Product Type: | CRISPRa Kits                                                                                                                                                                                                                                                                                            |
| Description:  | Sqor CRISPRa kit - CRISPR gene activation of mouse sulfide quinone oxidoreductase                                                                                                                                                                                                                       |
| Symbol:       | Sqor                                                                                                                                                                                                                                                                                                    |
| Components:   | GA206330G1, Sqor gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074)<br>GA206330G2, Sqor gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074)<br>GA206330G3, Sqor gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074)<br>1 CRISPRa-Enhancer vector, SKU GE100056<br>1 CRISPRa scramble vector, SKU GE100077 |
| UniProt ID:   | <u><a href="#">Q9R112</a></u>                                                                                                                                                                                                                                                                           |
| Synonyms:     | 0610039J17Rik; 4930557M22Rik; SQOR                                                                                                                                                                                                                                                                      |
| Format:       | 3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)                                                                                                                                                                                                                                 |

#### Product images:



[View online »](#)