

Product datasheet for **GA206363**

Mouse Tpm3 activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	Tpm3 CRISPRa kit - CRISPR gene activation of mouse tropomyosin 3, gamma
Symbol:	Tpm3
UniProt ID:	P21107
Synonyms:	gamma-TM; hTM30nm; hTMnm; Tm5NM; TM30nm; TMnm; Tpm-5; Tpm5; Trop-5
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online >](#)