

Product datasheet for **GA206378**

Mouse Myl10 activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	Myl10 CRISPRa kit - CRISPR gene activation of mouse myosin, light chain 10, regulatory
Symbol:	Myl10
Components:	GA206378G1, Myl10 gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074) GA206378G2, Myl10 gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074) GA206378G3, Myl10 gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074) 1 CRISPRa-Enhancer vector, SKU GE100056 1 CRISPRa scramble vector, SKU GE100077
UniProt ID:	<u>Q62082</u>
Synonyms:	1700027I08Rik; Mylc2pl; PLRLC
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online »](#)