

Product datasheet for **GA206384**

Mouse **Ms4a4b** activation kit by **CRISPRa**

Product data:

Product Type:	CRISPRa Kits
Description:	Ms4a4b CRISPRa kit – CRISPR gene activation of mouse membrane-spanning 4-domains, subfamily A, member 4B
Symbol:	Ms4a4b
Components:	GA206384G1, Ms4a4b gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GEI00074) GA206384G2, Ms4a4b gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GEI00074) GA206384G3, Ms4a4b gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GEI00074) 1 CRISPRa-Enhancer vector, SKU GEI00056 1 CRISPRa scramble vector, SKU GEI00077
Synonyms:	AI463180; Chandra; Ly116
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:

