

## Product datasheet for **GA206388**

### Mouse Sap30 activation kit by CRISPRa

#### Product data:

Product Type:	CRISPRa Kits
Description:	Sap30 CRISPRa kit - CRISPR gene activation of mouse sin3 associated polypeptide
Symbol:	Sap30
UniProt ID:	<u><a href="#">O88574</a></u>
Synonyms:	30kDa
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

#### Product images:



[View online >](#)