

Product datasheet for **GA206389**

Mouse Cenpk activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	Cenpk CRISPRa kit - CRISPR gene activation of mouse centromere protein K
Symbol:	Cenpk
Components:	GA206389G1, Cenpk gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074) GA206389G2, Cenpk gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074) GA206389G3, Cenpk gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074) 1 CRISPRa-Enhancer vector, SKU GE100056 1 CRISPRa scramble vector, SKU GE100077
UniProt ID:	Q9ESN5
Synonyms:	B130045K24Rik; BB232189; C530004N04Rik; Cenp-K; Solt; Solzt
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online »](#)