

Product datasheet for **GA206393**

Mouse Pgap6 activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	Tmem8 CRISPRa kit - CRISPR gene activation of mouse transmembrane protein 8
Symbol:	Pgap6
Components:	GA206393G1, Pgap6 gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074) GA206393G2, Pgap6 gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074) GA206393G3, Pgap6 gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074) 1 CRISPRa-Enhancer vector, SKU GE100056 1 CRISPRa scramble vector, SKU GE100077
UniProt ID:	Q9ESN3
Synonyms:	AW490096; M83; Tmem8a
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online »](#)