

Product datasheet for **GA206408**

Mouse Trp53inp1 activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	Trp53inp1 CRISPRa kit - CRISPR gene activation of mouse transformation related protein 53 inducible nuclear protein 1
Symbol:	Trp53inp1
Components:	GA206408G1, Trp53inp1 gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074) GA206408G2, Trp53inp1 gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074) GA206408G3, Trp53inp1 gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074) 1 CRISPRa-Enhancer vector, SKU GE100056 1 CRISPRa scramble vector, SKU GE100077
UniProt ID:	Q9QXE4
Synonyms:	2700057G22Rik; SIP; SIP18; SIP27; Stinp; Teap
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online »](#)