

## Product datasheet for **GA206410**

### Mouse *Kcnq4* activation kit by CRISPRa

#### Product data:

Product Type:	CRISPRa Kits
Description:	<i>Kcnq4</i> CRISPRa kit - CRISPR gene activation of mouse potassium voltage-gated channel, subfamily Q, member 4
Symbol:	<i>Kcnq4</i>
UniProt ID:	<a href="#">Q9JK97</a>
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

#### Product images:

