

## Product datasheet for **GA206413**

### Mouse *Taf8* activation kit by CRISPRa

#### Product data:

Product Type:	CRISPRa Kits
Description:	Taf8 CRISPRa kit - CRISPR gene activation of mouse TATA-box binding protein associated factor 8
Symbol:	Taf8
Components:	GA206413G1, Taf8 gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074) GA206413G2, Taf8 gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074) GA206413G3, Taf8 gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074) 1 CRISPRa-Enhancer vector, SKU GE100056 1 CRISPRa scramble vector, SKU GE100077
UniProt ID:	<a href="#">Q9EQH4</a>
Synonyms:	AW260255; Tbn
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

#### Product images:



[View online »](#)