

Product datasheet for **GA206415**

Mouse Impg1 activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	Impg1 CRISPRa kit - CRISPR gene activation of mouse interphotoreceptor matrix proteoglycan 1
Symbol:	Impg1
Components:	GA206415G1, Impg1 gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074) GA206415G2, Impg1 gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074) GA206415G3, Impg1 gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074) 1 CRISPRa-Enhancer vector, SKU GE100056 1 CRISPRa scramble vector, SKU GE100077
UniProt ID:	<u>Q8R1W8</u>
Synonyms:	A930015H12Rik; IMP150; SPACR
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online »](#)