

Product datasheet for **GA206434**

Mouse Gpr35 activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	Gpr35 CRISPRa kit - CRISPR gene activation of mouse G protein-coupled receptor 35
Symbol:	Gpr35
UniProt ID:	Q9ES90
Synonyms:	KPG_007
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online >](#)