


Product datasheet for **GA206438**

Mouse Sdf2l1 activation kit by CRISPRa

Product data:

| | |
|---------------|---|
| Product Type: | CRISPRa Kits |
| Description: | Sdf2l1 CRISPRa kit - CRISPR gene activation of mouse stromal cell-derived factor 2-like 1 |
| Symbol: | Sdf2l1 |
| Components: | GA206438G1, Sdf2l1 gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074) GA206438G2, Sdf2l1 gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074) GA206438G3, Sdf2l1 gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074) 1 CRISPRa-Enhancer vector, SKU GE100056 1 CRISPRa scramble vector, SKU GE100077 |
| UniProt ID: | <u>Q9ESP1</u> |
| Format: | 3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug) |

Product images:



[View online »](#)