

Product datasheet for **GA206463**

Mouse Ms4a4c activation kit by CRISPRa

Product data:

Product Type: CRISPRa Kits

Description: Ms4a4c CRISPRa kit - CRISPR gene activation of mouse membrane-spanning 4-domains, subfamily A, member 4C

Symbol: Ms4a4c

Components: GA206463G1, Ms4a4c gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GEI00074)
GA206463G2, Ms4a4c gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GEI00074)
GA206463G3, Ms4a4c gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GEI00074)
1 CRISPRa-Enhancer vector, SKU GEI00056
1 CRISPRa scramble vector, SKU GEI00077

Synonyms: 5830413L19Rik; AI449185; Ms4a9

Format: 3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:

