

Product datasheet for **GA206476**

Mouse Tspan4 activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	Tspan4 CRISPRa kit - CRISPR gene activation of mouse tetraspanin 4
Symbol:	Tspan4
UniProt ID:	Q9DCK3
Synonyms:	AI325509; AI746565; D130042I01Rik; NAG-2; Tm4sf7; Tspan-4
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online >](#)