

Product datasheet for **GA206482**

Mouse Mrps23 activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	Mrps23 CRISPRa kit - CRISPR gene activation of mouse mitochondrial ribosomal protein S23
Symbol:	Mrps23
Components:	GA206482G1, Mrps23 gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074) GA206482G2, Mrps23 gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074) GA206482G3, Mrps23 gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074) 1 CRISPRa-Enhancer vector, SKU GE100056 1 CRISPRa scramble vector, SKU GE100077
UniProt ID:	<u>Q8VE22</u>
Synonyms:	D11Bwg1153e; MRP-S23; Rpms23; S23mt
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online »](#)