


Product datasheet for **GA206489**

Mouse **Keg1** activation kit by **CRISPRa**

Product data:

Product Type:	CRISPRa Kits
Description:	Keg1 CRISPRa kit - CRISPR gene activation of mouse kidney expressed gene 1
Symbol:	Keg1
Components:	GA206489G1, Keg1 gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074) GA206489G2, Keg1 gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074) GA206489G3, Keg1 gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074) 1 CRISPRa-Enhancer vector, SKU GE100056 1 CRISPRa scramble vector, SKU GE100077
UniProt ID:	<u>Q9DCY0</u>
Synonyms:	0610008P16Rik; AI316519; GS4059
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online »](#)