

Product datasheet for **GA206496**

Mouse Lpin2 activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	Lpin2 CRISPRa kit - CRISPR gene activation of mouse lipin 2
Symbol:	Lpin2
Components:	GA206496G1, Lpin2 gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074) GA206496G2, Lpin2 gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074) GA206496G3, Lpin2 gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074) 1 CRISPRa-Enhancer vector, SKU GE100056 1 CRISPRa scramble vector, SKU GE100077
UniProt ID:	Q99PI5
Synonyms:	2610511G02Rik; AI481352; AW742896
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online »](#)