


Product datasheet for **GA206497**

Mouse Lpin3 activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	Lpin3 CRISPRa kit - CRISPR gene activation of mouse lipin 3
Symbol:	Lpin3
Components:	GA206497G1, Lpin3 gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074) GA206497G2, Lpin3 gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074) GA206497G3, Lpin3 gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074) 1 CRISPRa-Enhancer vector, SKU GE100056 1 CRISPRa scramble vector, SKU GE100077
UniProt ID:	<u>Q99PI4</u>
Synonyms:	9130206L11Rik; AA438110; AV236139; mKIAA4023
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online »](#)