

Product datasheet for **GA206507**

Mouse Lpar3 activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	Lpar3 CRISPRa kit - CRISPR gene activation of mouse lysophosphatidic acid receptor 3
Symbol:	Lpar3
UniProt ID:	<u>Q9EQ31</u>
Synonyms:	Edg7; lpA3
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online >](#)