

Product datasheet for **GA206539**

Mouse Tmem254a activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	Tmem254a CRISPRa kit - CRISPR gene activation of mouse transmembrane protein 254a
Symbol:	Tmem254a
Components:	GA206539G1, Tmem254a gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074) GA206539G2, Tmem254a gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074) GA206539G3, Tmem254a gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074) 1 CRISPRa-Enhancer vector, SKU GE100056 1 CRISPRa scramble vector, SKU GE100077
UniProt ID:	<u>P0DN89</u>
Synonyms:	0610008K04Rik; 0610010A22Rik; AA960404; D14Ert449e; Tmem254
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:

