

Product datasheet for **GA206554**

Mouse Snupn activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	Snupn CRISPRa kit - CRISPR gene activation of mouse snurportin 1
Symbol:	Snupn
Components:	GA206554G1, Snupn gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074) GA206554G2, Snupn gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074) GA206554G3, Snupn gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074) 1 CRISPRa-Enhancer vector, SKU GE100056 1 CRISPRa scramble vector, SKU GE100077
UniProt ID:	Q80W37
Synonyms:	0610031A09Rik; Rnut1; Snupn1
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online »](#)