

Product datasheet for **GA206561**

Mouse Tsen34 activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	Tsen34 CRISPRa kit - CRISPR gene activation of mouse tRNA splicing endonuclease subunit 34
Symbol:	Tsen34
Components:	GA206561G1, Tsen34 gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074) GA206561G2, Tsen34 gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074) GA206561G3, Tsen34 gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074) 1 CRISPRa-Enhancer vector, SKU GE100056 1 CRISPRa scramble vector, SKU GE100077
UniProt ID:	<u>Q8BMZ5</u>
Synonyms:	0610027F08Rik; Leng5
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online »](#)