


Product datasheet for **GA206571**

Mouse Lamtor4 activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	Lamtor4 CRISPRa kit - CRISPR gene activation of mouse late endosomal/lysosomal adaptor, MAPK and MTOR activator 4
Symbol:	Lamtor4
Components:	GA206571G1, Lamtor4 gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074) GA206571G2, Lamtor4 gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074) GA206571G3, Lamtor4 gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074) 1 CRISPRa-Enhancer vector, SKU GE100056 1 CRISPRa scramble vector, SKU GE100077
UniProt ID:	Q8CF66
Synonyms:	0910001L09Rik; AV006840
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online »](#)