

Product datasheet for **GA206580**

Mouse Tmed3 activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	Tmed3 CRISPRa kit - CRISPR gene activation of mouse transmembrane p24 trafficking protein 3
Symbol:	Tmed3
UniProt ID:	<u>Q78IS1</u>
Synonyms:	1200002G13Rik; AW546672; P24b
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online >](#)