


Product datasheet for **GA206595**

Mouse Tipin activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	Tipin CRISPRa kit - CRISPR gene activation of mouse timeless interacting protein
Symbol:	Tipin
Components:	GA206595G1, Tipin gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074) GA206595G2, Tipin gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074) GA206595G3, Tipin gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074) 1 CRISPRa-Enhancer vector, SKU GE100056 1 CRISPRa scramble vector, SKU GE100077
UniProt ID:	<u>Q91WA1</u>
Synonyms:	1110005A05Rik; 1110018P21Rik
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online »](#)