

Product datasheet for **GA206641**

Mouse Inip activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	Inip CRISPRa kit - CRISPR gene activation of mouse INTS3 and NABP interacting protein
Symbol:	Inip
Components:	GA206641G1, Inip gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074) GA206641G2, Inip gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074) GA206641G3, Inip gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074) 1 CRISPRa-Enhancer vector, SKU GE100056 1 CRISPRa scramble vector, SKU GE100077
UniProt ID:	<u>Q3TXT3</u>
Synonyms:	1110054O05Rik; 2610312O17Rik; AA399876; Ssbip1
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online »](#)