

Product datasheet for **GA206666**

Mouse Hsbp111 activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	Hsbp111 CRISPRa kit - CRISPR gene activation of mouse heat shock factor binding protein 1-like 1
Symbol:	Hsbp111
Components:	GA206666G1, Hsbp111 gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074) GA206666G2, Hsbp111 gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074) GA206666G3, Hsbp111 gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074) 1 CRISPRa-Enhancer vector, SKU GE100056 1 CRISPRa scramble vector, SKU GE100077
UniProt ID:	<u>B2RXB2</u>
Synonyms:	1810005K13Rik
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online »](#)