

Product datasheet for **GA206670**

Mouse Camk2n1 activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	Camk2n1 CRISPRa kit - CRISPR gene activation of mouse calcium/calmodulin-dependent protein kinase II inhibitor 1
Symbol:	Camk2n1
Components:	GA206670G1, Camk2n1 gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074) GA206670G2, Camk2n1 gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074) GA206670G3, Camk2n1 gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074) 1 CRISPRa-Enhancer vector, SKU GE100056 1 CRISPRa scramble vector, SKU GE100077
UniProt ID:	<u>Q6QWF9</u>
Synonyms:	1810006K23Rik; AI848599; AI891706
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online »](#)