

Product datasheet for **GA206676**

Mouse Pigyl activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	Pigyl CRISPRa kit - CRISPR gene activation of mouse phosphatidylinositol glycan anchor biosynthesis, class Y-like
Symbol:	Pigyl
Components:	GA206676G1, Pigyl gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074) GA206676G2, Pigyl gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074) GA206676G3, Pigyl gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074) 1 CRISPRa-Enhancer vector, SKU GE100056 1 CRISPRa scramble vector, SKU GE100077
UniProt ID:	<u>P0C1P0</u>
Synonyms:	1810008A14Rik; PIG-Y; Pigy
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online »](#)