

Product datasheet for **GA206677**

Mouse Tmed6 activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	Tmed6 CRISPRa kit - CRISPR gene activation of mouse transmembrane p24 trafficking protein 6
Symbol:	Tmed6
Components:	GA206677G1, Tmed6 gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074) GA206677G2, Tmed6 gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074) GA206677G3, Tmed6 gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074) 1 CRISPRa-Enhancer vector, SKU GE100056 1 CRISPRa scramble vector, SKU GE100077
UniProt ID:	Q9CQG0
Synonyms:	1810015P03Rik; 1810018I24Rik
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online »](#)