

Product datasheet for **GA206687**

Mouse Tma16 activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	Tma16 CRISPRa kit - CRISPR gene activation of mouse translation machinery associated 16
Symbol:	Tma16
Components:	GA206687G1, Tma16 gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074) GA206687G2, Tma16 gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074) GA206687G3, Tma16 gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074) 1 CRISPRa-Enhancer vector, SKU GE100056 1 CRISPRa scramble vector, SKU GE100077
UniProt ID:	<u>Q9CR02</u>
Synonyms:	1810029B16Rik
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online »](#)