

Product datasheet for **GA206717**

Mouse Psenen activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	Psenen CRISPRa kit - CRISPR gene activation of mouse presenilin enhancer gamma secretase subunit
Symbol:	Psenen
UniProt ID:	Q9CQR7
Synonyms:	1700023M09Rik; Pen-2
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online >](#)