

## Product datasheet for **GA206771**

### Mouse **Pcp4l1** activation kit by **CRISPRa**

#### Product data:

|               |                                                                                                                                                                                                                                                                                                               |
|---------------|---------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| Product Type: | CRISPRa Kits                                                                                                                                                                                                                                                                                                  |
| Description:  | Pcp4l1 CRISPRa kit - CRISPR gene activation of mouse Purkinje cell protein 4-like 1                                                                                                                                                                                                                           |
| Symbol:       | Pcp4l1                                                                                                                                                                                                                                                                                                        |
| Components:   | GA206771G1, Pcp4l1 gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074)<br>GA206771G2, Pcp4l1 gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074)<br>GA206771G3, Pcp4l1 gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074)<br>1 CRISPRa-Enhancer vector, SKU GE100056<br>1 CRISPRa scramble vector, SKU GE100077 |
| UniProt ID:   | <a href="#">Q6W8Q3</a>                                                                                                                                                                                                                                                                                        |
| Synonyms:     | 1500036F01Rik                                                                                                                                                                                                                                                                                                 |
| Format:       | 3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)                                                                                                                                                                                                                                       |

#### Product images:



[View online »](#)