

Product datasheet for **GA206818**

Mouse Psmg3 activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	Psmg3 CRISPRa kit - CRISPR gene activation of mouse proteasome (prosome, macropain) assembly chaperone 3
Symbol:	Psmg3
Components:	GA206818G1, Psmg3 gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074) GA206818G2, Psmg3 gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074) GA206818G3, Psmg3 gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074) 1 CRISPRa-Enhancer vector, SKU GE100056 1 CRISPRa scramble vector, SKU GE100077
UniProt ID:	<u>Q9CZH3</u>
Synonyms:	1810042K04Rik; 4930403H09Rik; AI303239
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online »](#)