

Product datasheet for **GA206865**

Mouse Ms4a4d activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	Ms4a4d CRISPRa kit – CRISPR gene activation of mouse membrane-spanning 4-domains, subfamily A, member 4D
Symbol:	Ms4a4d
Components:	GA206865G1, Ms4a4d gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GEI00074) GA206865G2, Ms4a4d gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GEI00074) GA206865G3, Ms4a4d gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GEI00074) 1 CRISPRa-Enhancer vector, SKU GEI00056 1 CRISPRa scramble vector, SKU GEI00077
UniProt ID:	<u>Q99N05</u>
Synonyms:	2210417J23Rik
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:

