

Product datasheet for **GA206877**

Mouse Thg1l activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	Thg1l CRISPRa kit - CRISPR gene activation of mouse tRNA-histidine guanylyltransferase 1-like (<i>S. cerevisiae</i>)
Symbol:	Thg1l
Components:	GA206877G1, Thg1l gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074) GA206877G2, Thg1l gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074) GA206877G3, Thg1l gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074) 1 CRISPRa-Enhancer vector, SKU GE100056 1 CRISPRa scramble vector, SKU GE100077
UniProt ID:	<u>Q9CY52</u>
Synonyms:	1700121M19Rik; 5730409G07Rik; AA387658; RP23-298M7.6
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online »](#)