

## Product datasheet for **GA206881**

### Mouse Tsen15 activation kit by CRISPRa

#### Product data:

Product Type:	CRISPRa Kits
Description:	Tsen15 CRISPRa kit - CRISPR gene activation of mouse tRNA splicing endonuclease subunit 15
Symbol:	Tsen15
UniProt ID:	<a href="#">Q8R3W5</a>
Synonyms:	5730449L18Rik; AL023077; Sen15
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

#### Product images:



[View online >](#)