

Product datasheet for **GA206896**

Mouse Msantd3 activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	Msantd3 CRISPRa kit - CRISPR gene activation of mouse Myb/SANT-like DNA-binding domain containing 3
Symbol:	Msantd3
Components:	GA206896G1, Msantd3 gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074) GA206896G2, Msantd3 gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074) GA206896G3, Msantd3 gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074) 1 CRISPRa-Enhancer vector, SKU GE100056 1 CRISPRa scramble vector, SKU GE100077
UniProt ID:	<u>Q9CR78</u>
Synonyms:	2410046L14Rik; 2810429P03Rik; 5730528L13Rik; Mysantd3
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online »](#)