


Product datasheet for **GA206901**

Mouse Tmed7 activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	Tmed7 CRISPRa kit - CRISPR gene activation of mouse transmembrane p24 trafficking protein 7
Symbol:	Tmed7
Components:	GA206901G1, Tmed7 gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074) GA206901G2, Tmed7 gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074) GA206901G3, Tmed7 gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074) 1 CRISPRa-Enhancer vector, SKU GE100056 1 CRISPRa scramble vector, SKU GE100077
Synonyms:	3930401E15Rik; 5830493P14Rik; AU044611; AW545765; TRAM
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online »](#)