

Product datasheet for **GA206930**

Mouse *Kcng4* activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	<i>Kcng4</i> CRISPRa kit - CRISPR gene activation of mouse potassium voltage-gated channel, subfamily G, member 4
Symbol:	<i>Kcng4</i>
UniProt ID:	Q80XM3
Synonyms:	4921535I01Rik; AW049024; KV6.3; KV6.4
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online >](#)