

Product datasheet for GA206977

Mouse Ttc1 activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	Ttc1 CRISPRa kit - CRISPR gene activation of mouse tetratricopeptide repeat domain 1
Symbol:	Ttc1
Components:	GA206977G1, Ttc1 gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074) GA206977G2, Ttc1 gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074) GA206977G3, Ttc1 gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074) 1 CRISPRa-Enhancer vector, SKU GE100056 1 CRISPRa scramble vector, SKU GE100077
UniProt ID:	Q91Z38
Synonyms:	4833412C19Rik; C79328; TPR1
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:

